Factor 10a (Xa) Inhibition, Low Molecular Weight Heparin (LMWH)

Order Name: LMWHEPARIN

Test Number: 1506175
Revision Date: 08/26/2014

TEST NAME			METHO	DOLOGY	LOINC CODE
Factor 10a (Xa) Inhibition, Low Molecular Weight Heparin (LMWH)			Chromogenic		32684-3
SPECIMEN REQUIREM	MENTS				
Specimen	Specimen Volume (min)	Specimen Type		Specimen Container	Transport Environment
Preferred	2.7 mL	Whole Blood		Sodium Citrate 3.2% (Blue Top)	Room Temperature
Alternate 1	1.5 mL	Double Spun Plasma		Sterile, Capped Plastic Tube	Frozen
Instructions	Please indicate anticoagulant therapy. Each 2.7mL Sodium Citrate 3.2% (Blue Top) tube must be filled to the proper level, no hemolysis. Improperly filled tubes can give erroneous results. Specimen should be collected at least 4 hours after subcutaneous injection. Whole blood must be transported to SJMC lab within 1 hour of collection. If testing cannot be started within 1 hour of collection the specimen must be double spun then 1.5mL plasma aliquot from each tube into individual plastic aliquot tubes and freeze. Do not pool aliquots together!				

GENERAL INFORMATION		
Testing Schedule	Daily	
Expected TAT	1 Day	
Clinical Use	The assesment of the anti-Xa effect of heparin by this method gives a high sensitivity in detecting Low Molecular Weight Heparin levels.	
CPT Code(s)	85520	
Lab Section	Coagulation	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025