

## Inhibin A

Order Name: **INHIBIN A**  
Test Number: 3622375  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Inhibin A	Immunochemiluminometric (ICMA)	23883-2

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>0.3 mL (0.2)</b>	<b>Serum</b>	<b>Clot Activator SST</b>	<b>Frozen</b>
Alternate 1	<b>0.3 mL (0.2)</b>	<b>Serum</b>	<b>Clot Activator (Red Top, No-Gel)</b>	<b>Frozen</b>
Instructions	<p><b>Notes:</b> 0.2 mL (Note: This volume Does NOT allow for repeat testing.) <b>Specimen Type:</b> Red-top tube or gel-barrier tube <b>Specimen Storage:</b> <b>Freeze</b> <b>Specimen Collection:</b> Transfer the serum into a LabCorp PP transpak frozen purple tube with screw cap (LabCorp ID 23751). The specimen should be FROZEN immediately and maintained frozen until tested. To avoid delays in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED. <b>Special Instructions:</b> This test may exhibit interference when sample is collected from a person who is consuming a supplement with a high dose of biotin (also termed as vitamin B7 or B8, vitamin H, or coenzyme R). It is recommended to ask all patients who may be indicated for <b>Specimen Stability:</b> Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p>			

### GENERAL INFORMATION

Expected TAT	4 - 6 days
Performing Labcorp Test Code	146803
Notes	Labcorp Test Code: 146803
CPT Code(s)	86336
Lab Section	Reference Lab