

## DNase-B Antibody

Order Name: **DNASE B AB**  
Test Number: 5500210  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
DNase-B Antibody	Immunonephelometry	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>1 mL (0.5)</b>	<b>Serum</b>	<b>Clot Activator SST</b>	<b>Room Temperature</b>
Alternate 1	<b>1 mL (0.5)</b>	<b>Serum</b>	<b>No Additive Clot (Red Top, No-Gel, Plastic)</b>	<b>Room Temperature</b>

**Instructions**  
**Specimen Type:** Red-top tube or gel-barrier tube  
**Specimen Storage:** Maintain specimen at room temperature.  
**Specimen Collection:** Not Available  
**Specimen Stability:** Ambient: 9 days, Refrigerated : 14 days, Frozen: 14 days

### GENERAL INFORMATION

Expected TAT	3 - 5 days
Clinical Use	DNase-B Antibody is useful in patients with group A Streptococcal infection. DNase-B Antibody may persist for as long as three months. Comparison of titers of acute and convalescent specimens is useful for diagnosis of group A streptococcal infection.
Performing Labcorp Test Code	096289
Notes	Labcorp Test Code: 096289
CPT Code(s)	86215
Lab Section	Reference Lab