

## Hepatitis C Analyzer Reflex to Genotype

Order Name: **HEP C AN+G**  
Test Number: 5590857  
Revision Date: 05/28/2021

TEST NAME	METHODOLOGY	LOINC CODE
Hepatitis C Analyzer Reflex to Genotype	Chemiluminescence Assay	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<a href="#">See Instructions</a>	<b>Serum and Plasma</b>	<a href="#">See Instructions</a>	<a href="#">See Instructions</a>
Instructions	Two Specimens are Required for this Analyzer.			
<p><b>Specimen 1</b> - for the Hepatitis C Antibody  <b>4mL(2mL) Serum from Clot Activator SST, Keep Refrigerated</b> or Plasma from Lithium Heparin PST (Light Green)          Stability: Ambient 3 days, refrigerated 7 days, Frozen &gt;7 days.</p> <p><b>Specimen 2</b> - for possible HCV Quantitative Viral Load by PCR  <b>Preferred Specimen is 2.5 mL (1.5 mL) Plasma or Serum from EDTA (Lavender Top) or Clot Activator SST.</b>  <b>Centrifuge specimen and separate plasma or serum from cells within 2 hours of collection.</b>          Alternately specimens may be held at 15 to 30°C for up to 6 hours prior to separation and centrifugation.          After centrifugation, serum or plasma may be held for up to 72 hours at 2-8°C.          Serum or plasma specimens may be stored frozen at -20°C or lower for up to 60 days. Avoid multiple freeze/thaws.</p>				

### GENERAL INFORMATION

Testing Schedule	Mon-Fri
Expected TAT	1-2 Days (Assay Dependant)
Clinical Use	Positive Hepatitis C Antibody will order the following test for viral load and genotype determination. <b>Hepatitis C RNA Quantitative PCR with reflex to Genotyping [5594677] HCV Q+GENO</b>
Notes	Assay Schedule: HCV Quantitative: Dayshift- Set up on Monday, Wednesday and Friday reported same day. HCV Genotype: Dayshift- Set up on Thursday and reported on Friday.
CPT Code(s)	86803 (possible 87522 for HCV Quantitative Viral Load by PCR ) (possible 87902 for HCV Genotyping )
Lab Section	Chemistry