

HCV GenoSure(R) NS3/4A

Order Name: **Hep C NS3/4A**

Test Number: 5194390

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
HCV GenoSure(R) NS3/4A	Polymerase Chain Reaction/Sequencing	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1 mL)	Plasma	EDTA (Lavender Top)	Frozen
Instructions	<p><b>Specimen Type:</b> Ship in plasma preparation tube (PPT(TM)) OR screw-capped polypropylene frozen transport tube.</p> <p><b>Specimen Storage:</b> FREEZE</p> <p><b>Specimen Collection:</b> Collect specimen in one PPT(TM) tube, gel-barrier tube, or lavender-top (EDTA) tube. Do NOT use green-top (heparin) tubes. Centrifuge specimen within six hours of collection, remove serum/plasma and transfer specimen to polypropylene screw-capped tubes, and FREEZE. Specimens drawn into PPT(TM) tubes need to be centrifuged within six hours and frozen. Ship FROZEN. To avoid delays in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.</p> <p><b>Special Instructions:</b> Please indicate the patient's HCV Genotype (1a or 1b) on the test request form. Testing referred to Monogram Biosciences L9#C6000</p> <p><b>Specimen Stability:</b> Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p>			

GENERAL INFORMATION	
Expected TAT	10 - 14 days
Performing Labcorp Test Code	550540
Notes	Polymerase chain reaction (PCR) amplification and DNA next generation sequencing Labcorp Test Code: 550540
CPT Code(s)	87902, 87900
Lab Section	Reference Lab