

Microalbumin Urine Timed

Order Name: **ALBUM24U**
Test Number: 2022250
Revision Date: 05/28/2021

| TEST NAME | METHODOLOGY | LOINC CODE |
|--|----------------------------------|------------|
| Total Urine Volume | | 3167-4 |
| Creatinine, Urine, mg/dL | Kinetic Alkaline Picrate (Jaffe) | 20624-3 |
| Microalbum/Creat Ratio 24H UR | | 14958-3 |
| Microalbumin Excretion | | 58448-2 |
| Microalbuminuria 24 Hour Urine mg/24hr | | 14956-7 |
| Microalbuminuria 24 Hour Urine mg/l | Immunoturbidimetry | 30003-8 |
| Hours Collected - Prompt | | 30211-7 |
| Volume 24 Urine Prompt | | |

SPECIMEN REQUIREMENTS

| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
|--------------|---|----------------|-------------------------|-----------------------|
| Preferred | 1 mL (0.5) | Urine, 24-hour | 24 hour Urine Container | Refrigerated |
| Instructions | 24 hour urine collection. No preservative. (6N HCl preservative is acceptable if collecting with other tests.) Record number of hours and volume in ml on the specimen container. Keep refrigerated. Specimen stability: Ambient 24hrs, Refrigerated 6 days Unacceptable Specimens - Grossly turbid or bloody urine | | | |

GENERAL INFORMATION

| | |
|------------------|--|
| Testing Schedule | Mon - Fri |
| Expected TAT | 1-2 days |
| Clinical Use | Useful for diabetic patients to assess the potential for early onset of nephropathy. |
| CPT Code(s) | 82043 |
| Lab Section | Chemistry |