

Metanephrines, Plasma (Free)

Order Name: **METANEPH P**
 Test Number: 3805327
 Revision Date: 02/16/2021

TEST NAME	METHODOLOGY	LOINC CODE
Metanephrines, Plasma (Free)	Liquid Chromatography/Tandem Mass Spectrometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.6 mL (0.4 mL)	Plasma	EDTA (Lavender Top)	Refrigerated
Instructions	<p>Specimen Type: Lavender-top (EDTA) tube</p> <p>Specimen Storage: Sample must be refrigerated due to limited stability at room temperature. Freeze if unable to maintain refrigerated.</p> <p>Specimen Collection: Draw blood in chilled lavender-top (EDTA) tube. Invert to mix with preservatives. Centrifuge and transfer the plasma to a labeled plastic transport tube. Refrigerate or freeze separated plasma immediately. IMPORTANT Note: The patient should be in a fully recumbent position (lying down) for at least 20 minutes before and during sample collection.</p> <p>Special Instructions: For those centers where sampling blood in the fully recumbent, supine position is not possible testing for 24-hour urine metanephrines(LabCorp test number test 004234) may be preferable.</p> <p>Specimen Stability: Ambient: Unstable, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	3-6 Days
Performing Labcorp Test Code	121806
Notes	Labcorp Test Code: 121806
CPT Code(s)	83835
Lab Section	Reference Lab