

## Candida Antibody

Order Name: **CANDIDA AB**  
Test Number: 5569300  
Revision Date: 07/10/2025

| TEST NAME        | METHODOLOGY                  | LOINC CODE |
|------------------|------------------------------|------------|
| Candida Antibody | Double immunodiffusion (DID) | 5074-0     |

### SPECIMEN REQUIREMENTS

|              |   |               |   |                         |
|--------------|---|---------------|---|-------------------------|
| Specimen     | Specimen Volume (min)   | Specimen Type | Specimen Container                        | Transport Environment   |
| Preferred    | <b>1 mL (0.5mL)</b>   | <b>Serum</b>  | <b>Clot Activator (SST or Red No-Gel)</b> | <b>Room Temperature</b> |
| Instructions | <b>Specimen:</b> 1mL (0.5mL Note: This volume does not allow for repeat testing.) Serum from Red-top tube or gel-barrier tube<br><b>Cause for Rejection:</b> Lipemia; gross bacterial contamination |               |   |                         |

### GENERAL INFORMATION

|                              |   |
|------------------------------|---|
| Expected TAT                 | 2-8 Days  |
| Clinical Use                 | The literature currently has conflicting data on the usefulness of serologic testing for candida antibodies in patients with candidiasis. |
| Performing Labcorp Test Code | 096719  |
| CPT Code(s)                  | 86628   |
| Lab Section                  | Reference Lab   |