

Urea Nitrogen Urine Random

Order Name: **UUN R U**

Test Number: 3000270

Revision Date: 06/04/2021

| TEST NAME | METHODOLOGY | LOINC CODE |
|----------------------------|-------------|------------|
| Urea Nitrogen Urine Random | Urease/GLDH | 3095-7 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|-------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.5) | Urine, Random | Sterile Urine container | Refrigerated |
| Instructions | Random urine collection. No preservative. Keep refrigerated. Specimen stability: Ambient 2 days, Refrigerated 7 days, Frozen 1 month | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Testing Schedule | Daily |
| Expected TAT | 1-2 days |
| Clinical Use | Useful as an aid in evaluating renal function. |
| CPT Code(s) | 84540 |
| Lab Section | Chemistry |