

## Follicle Stimulating Hormone - FSH, Serum

Order Name: **FSH**  
Test Number: 3601200  
Revision Date: 05/20/2021

TEST NAME	METHODOLOGY	LOINC CODE
Follicle Stimulating Hormone - FSH, Serum	Chemiluminescence Assay	15067-2

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>1 mL (0.5)</b>	<b>Serum</b>	<b>Clot Activator SST</b>	<b>Refrigerated</b>
Alternate 1	<b>1 mL (0.5)</b>	<b>Plasma</b>	<b>Lithium Heparin PST (Light Green Top)</b>	<b>Refrigerated</b>
Instructions	Remove Serum/Plasma from cells/gel if testing delayed >24hrs. Stability: Ambient 1 day. Refrigerated 7 days. Frozen 12 months.			

### GENERAL INFORMATION

Testing Schedule	Mon - Fri
Expected TAT	1-2 days
Clinical Use	Useful for evaluation of menstrual irregularities, work-up of patients with suspected hypogonadism, prediction of ovulation, evaluation of infertility and the diagnosis of pituitary disorders.
CPT Code(s)	83001
Lab Section	Chemistry