

Acetaminophen Quantitative

Order Name: **ACETAMIN**

Test Number: 4000050

Revision Date: 04/26/2021

| TEST NAME | METHODOLOGY | LOINC CODE |
|----------------------------|-------------|------------|
| Acetaminophen Quantitative | Enzymatic | 3298-7 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|---------------------------------------|-----------------------|
| Specimen | Specimen Volume (mL) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.5) | Plasma | Lithium Heparin PST (Light Green Top) | Refrigerated |
| Alternate 1 | 1 mL (0.5) | Serum | Clot Activator SST | Refrigerated |
| Instructions | Recommended collection time is four hours after an oral dose. Specimen stability: Ambient 24 hrs, Refrigerated 14 days, Frozen 45 days | | | |

| GENERAL INFORMATION | |
|---------------------|---|
| Testing Schedule | Daily |
| Expected TAT | 1-2 days |
| Clinical Use | Useful for monitoring toxicity in overdose cases. |
| CPT Code(s) | 80329 |
| Lab Section | Chemistry |