

Apolipoprotein B

Order Name: **APO B**
Test Number: 2015050
Revision Date: 04/26/2021

| TEST NAME | METHODOLOGY | LOINC CODE |
|------------------|---------------|------------|
| Apolipoprotein B | Turbidometric | 1871-3 |

SPECIMEN REQUIREMENTS

| | | | | |
|--------------|---|---------------|--|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.5) | Serum | Clot Activator SST | Refrigerated |
| Alternate 1 | 1 mL (0.5) | Plasma | Lithium Heparin PST (Light Green Top) | Refrigerated |
| Instructions | <p>Patient Preparation: Freshly drawn fasting specimen. Separate serum/plasma from cells ASAP or within 2 hours of collection Unacceptable Conditions: Hemolyzed specimens. Stability After separation from cells: Ambient 8 hours, Refrigerated 7 days, Frozen 2 months</p> | | | |

GENERAL INFORMATION

| | |
|------------------|-----------|
| Testing Schedule | Mon - Fri |
| Expected TAT | 2-3 Days |
| CPT Code(s) | 82172 |
| Lab Section | Chemistry |