

Beta-2, Transferrin BF

Order Name: **BETA-2 TRA**

Test Number: 3656675

Revision Date: 10/30/2025

TEST NAME	METHODOLOGY	LOINC CODE
Beta-2, Transferrin BF	Immunofixation Electrophoresis	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL	Body Fluid	Sterile Screwtop Container	Refrigerated
Instructions	Aural or nasal fluid in sterile container; Do Not Freeze. Do Not send CSF specimens! Notes: 1 mL (Note: This volume Does NOT allow for repeat testing). Specimen Type: Sterile container Specimen Storage: Refrigerated Special Instructions: Testing referred to ARUP TC 50047 Specimen Stability: Ambient: 2 weeks; Refrigerated: 2 weeks; Frozen: 2 weeks			

GENERAL INFORMATION	
Expected TAT	4 - 8 days
Clinical Use	This test is to detect CSF in body fluids such as sinus or ear drainage.
Performing Labcorp Test Code	829030
CPT Code(s)	86335
Lab Section	Reference Lab