

KOH Prep for Fungus

Order Name: **C KOH PR**

Test Number: 6300100

Revision Date: 02/21/2020

| TEST NAME | METHODOLOGY | LOINC CODE |
|---------------------|-------------|------------|
| KOH Prep for Fungus | Microscopy | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|------------------|--------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | See Instructions | See Instructions | See Instructions | Room Temperature |
| Instructions | Submit tissue, skin scrapings, or nail pieces in a sterile screw-top cup. For vaginal microscopic exam, send vaginal discharge in a saline wet prep tube. | | | |

| GENERAL INFORMATION | |
|---------------------|---|
| Testing Schedule | Daily |
| Expected TAT | 1 Day |
| Clinical Use | Detection of fungal elements in tissue samples; detection of yeast/fungal elements in vaginal samples |
| CPT Code(s) | 87220 |
| Lab Section | Microbiology |