

Barbiturate Screen, Urine

Order Name: **UBARBS**  
Revision Date: 01/02/2020

TEST NAME		METHODOLOGY		LOINC CODE
Barbiturate Screen, Urine				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	7 mL (1 mL)	Urine, Random	Sterile Container, No Preservative	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided			
Methodology	Homogeneous enzyme immunoassay			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	3 Hours
Stat TAT	45 Minutes
Notes	Specimen Stability: Room Temp - 7 Days Refrigerated - 7 Days Frozen - Indefinite
CPT Code(s)	G0431
Lab Section	NRLS-Core Chemistry