

D-Dimer

Order Name: **DD**
Revision Date: 11/13/2019

TEST NAME		METHODOLOGY		LOINC CODE
D-Dimer				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.8 mL (1.8 mL)	Plasma	3.2% Sodium Citrate tube	Room Temperature
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			

GENERAL INFORMATION	
Testing Schedule	Monday through Sunday
Expected TAT	3 Hours
Stat TAT	1 Hour
Notes	Specimen Stability: Ambient - 4 Hours If testing cannot be performed within 4 hours (2 hours if opened), freeze an aliquot of the plasma.
CPT Code(s)	85379
Lab Section	NRLS-Core Hematopathology