

Gentamicin, Peak

Order Name: **GENP**
Revision Date: 11/19/2019

| TEST NAME | | METHODOLOGY | | LOINC CODE |
|-----------------------|---|---------------|--------------------|-----------------------|
| Gentamicin, Peak | | | | |
| SPECIMEN REQUIREMENTS | | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1.0 mL (0.5 mL) | Serum | Serum Gel Tube | Refrigerated |
| Instructions | Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered. | | | |
| Reference Range | An interpretive report will be provided. | | | |
| Methodology | Homogeneous Enzyme Immunoassay | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Testing Schedule | Monday through Sunday |
| Expected TAT | 3 Hours |
| Stat TAT | 45 Minutes |
| Notes | Specimen Stability: Ambient: 8 Hours Refrigerated: 6 Weeks Frozen: 6 Months |
| CPT Code(s) | 80170 |
| Lab Section | NRLS-Core Chemistry |