

Immunoglobulin A

Order Name: **IGA**
Revision Date: 01/12/2013

TEST NAME	METHODOLOGY	LOINC CODE
Immunoglobulin A		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Instructions	Draw blood in a serum gel tube(s). Spin down and send 1 mL (minimum volume: 0.5 mL) of serum refrigerated. Note: Label specimen with patient name (first and last) or other unique identifier, patient hospital identification number or Social Security number, date and time of draw, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			

GENERAL INFORMATION	
Testing Schedule	Monday through Sunday
CPT Code(s)	82784
Lab Section	NRLS-Core Chemistry