

Opiate Screen, Urine

Order Name: **UOPIS**
Revision Date: 01/12/2013

TEST NAME		METHODOLOGY		LOINC CODE
Opiate Screen, Urine				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	7 mL (1 mL)	Urine, Random	Preservative-free, clean container	Refrigerated
Instructions	This test is for non-legal purposes for medical emergency only. All other requests are sent to St. Anthony Toxicology Laboratory. Submit 7 mL (minimum volume: 1 mL) from a random urine collection in a clean, dry, screw-capped container. No preservative. Send specimen refrigerated. Note: Label specimen with patient name (first and last) or other unique identifier, patient hospital identification number or Social Security number, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	See Interpretive Report			
Methodology	Enzyme Immunoassay			

GENERAL INFORMATION	
Testing Schedule	Monday through Sunday
Expected TAT	3 Hours
Stat TAT	45 Minutes
CPT Code(s)	80101
Lab Section	NRLS-Core Chemistry