

Rapid Plasma Reagin (RPR)

Order Name: **RPR**  
Revision Date: 01/12/2013

TEST NAME	METHODOLOGY	LOINC CODE
Rapid Plasma Reagin (RPR)		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Instructions	Draw blood in a serum gel tube(s). Spin down and send 1 mL (minimum volume: 0.2 mL) of serum refrigerated. Note: Label specimen with patient name (first and last) or other unique identifier, patient hospital identification number or Social Security number, date and time of draw, collector initials, and test(s) being ordered.			
Reference Range	Nonreactive			

GENERAL INFORMATION	
Testing Schedule	Monday through Friday (if received Monday-Thursday) TAT: 24 hours
CPT Code(s)	86592
Lab Section	NRLS-Microbiology