

Vaginosis Panel

Order Name: **VAGDNAP**  
Revision Date: 06/16/2014

TEST NAME		METHODOLOGY		LOINC CODE
Vaginosis Panel				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Reference Range	No Candida species, Gardnerella vaginalis , or Trichomonas vaginalis detected			
Methodology	DNA Probe Includes Candida species, Gardnerella vaginalis , and Trichomonas vaginalis .			

GENERAL INFORMATION	
Testing Schedule	Monday through Sunday
Clinical Use	<div>1. Obtain specimen from vagina site using BD Affirm VPIII Swab Collection Kit.</div> <div>2. Insert swab into vagina and rotate for adequate sampling.</div> <div>3. Place swab into BD Affirm VP111 Transport System.</div> <div>Note: Samples transported in VPIII Transport System are stable for up to 72 hours. The system has also been qualified for refrigerated transport.</div> <div>4. Label tube with patient's name (first and last) or other unique identifier, patient's hospital identification number or Social Security number, date and time of collection, collector's initials, and test(s) being ordered.</div> <div>5. Maintain sterility and forward promptly at ambient temperature only.</div> <div>Note: Specimen source and suspected species are required on request form for processing.</div>
CPT Code(s)	87480 - Candida species 87510 - Gardnerella vaginalis 87660 - Trichomonas vaginalis
Lab Section	NRLS-Microbiology