

WBC Count

Order Name: **WBC**  
Revision Date: 01/12/2013

TEST NAME		METHODOLOGY		LOINC CODE
WBC Count				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Instructions	Draw blood in a lavender-top (K2 EDTA) tube(s), and send 4 mL (minimum volume: 1 mL) of K2 EDTA whole blood refrigerated. Note: Label specimen with patient name (first and last) or other unique identifier, patient hospital identification number or Social Security number, date and time of draw, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			
GENERAL INFORMATION				
Testing Schedule	Monday through Sunday TAT: 1 hour STAT TAT: 30 minutes			
CPT Code(s)	85048			
Lab Section	NRLS-Core Hematopathology			