## Elm Tree Allergen

Order Name:ELM TREETest Number:5606795Revision Date:02/11/2013

| TEST NAME  |  |               | METHODOLOGY        | LOINC CODE               |
|--|--|---------------|--------------------|--------------------------|
| Elm Tree Allergen  |  |               | ImmunoCAP          | 15697-6                  |
| SPECIMEN REQUI   | REMENTS  |               |                    |                          |
| Specimen   | Specimen Volume (min)  | Specimen Type | Specimen Container | Transport Environment    |
| Preferred  | 0.4 mL (0.2 mL)  | Serum         | Clot Activator SST | Room Temperature         |
| Instructions   | Specimen Type: Red-top tube or gel-barrier tube,<br>Separate Serum from Cells into a screwtop transport container. |               |                    |                          |
| Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze |  |               |                    | e/thaw cycles Stable x3) |

| GENERAL INFORMATION |               |  |  |
|---------------------|---------------|--|--|
| Expected TAT        | 3-5 days      |  |  |
| Notes               |               |  |  |
| CPT Code(s)         | 86003         |  |  |
| Lab Section         | Reference Lab |  |  |