

Meconium Drug Screen (5 Panel)

Order Name: **DRUG F SCR**

Test Number: 4066250

Revision Date: 12/20/2022

| TEST NAME | METHODOLOGY | LOINC CODE |
|--------------------------------|------------------|------------|
| Meconium Drug Screen (5 Panel) | Immunoassay (IA) | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|------------------|--------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 5g (1g) | See Instructions | See Instructions | Refrigerated |
| Instructions | CHAIN OF CUSTODY TEST ! Specimen Type: Meconium Collection Kit (PeopleSoft#96663) OR clean, leakproof polypropylene container. Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen Collection: Specimens from different voidings may be pooled if necessary Special Instructions: Testing referred to MedTox Laboratories Inc MX#4034 | | | |

| GENERAL INFORMATION | |
|------------------------------|---|
| Expected TAT | 3-6 Days |
| Clinical Use | All fecal drug screens are analyzed by immunoassay. Each positive result is confirmed by GC/MS or LC/MS/MS to the limit of detection. Meconium Drug Screen Confirmations |
| Performing Labcorp Test Code | 809373 |
| Notes | Labcorp Test Code: 809373 |
| CPT Code(s) | 80307 |
| Lab Section | Reference Lab |