

Diuretics Screen, Urine

Order Name: **DIURET S U**

Test Number: 3632070

Revision Date: 01/04/2023

TEST NAME	METHODOLOGY	LOINC CODE
Diuretics Screen, Urine	Liquid Chromatography/Tandem Mass Spectrometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (1.2 mL)	Urine, Random	Sterile Screwtop Container	Refrigerated
Instructions	<p>Specimen Type: Plastic urine container without preservative</p> <p>Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p>Special Instructions: Testing referred to MedTox Laboratories Inc MX#413</p>			

GENERAL INFORMATION	
Expected TAT	10-14 Days
Performing Labcorp Test Code	808548
Notes	Labcorp Test Code: 808548
CPT Code(s)	80377
Lab Section	Reference Lab