

Azithromycin Level by HPLC, Serum or Plasma

Order Name: **Azithromycin Lvl**

Test Number: 6904757

Revision Date: 01/02/2018

| TEST NAME | METHODOLOGY | LOINC CODE |
|---|--|------------|
| Azithromycin Level by HPLC, Serum or Plasma | Quantitative HPLC/Tandem Mass Spectrometry | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|----------------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 2 mL | Serum | Clot Activator (Red Top, No-Gel) | Frozen |
| Instructions | <p>Specimen Type: Red-top tube</p> <p>Specimen Storage: FREEZE</p> <p>Specimen Collection: Serum should be separated from cells within two hours of venipuncture. Send serum in a plastic transport tube before freezing.</p> <p>Special Instructions: Testing referred to National Jewish Health NJHCO#PKAZM Was AZL NOTE: Drug dose and frequency, including time of last dose should be indicated on request form.</p> | | | |

| GENERAL INFORMATION | |
|------------------------------|---------------------------|
| Expected TAT | 11-17 days |
| Performing Labcorp Test Code | 857340 |
| Notes | Labcorp Test Code: 857340 |
| CPT Code(s) | 80299 |
| Lab Section | Reference Lab |