

Arbovirus Antibodies, IgG and IgM, Serum

Order Name: **Arbovirus Ab IgG/IgM**

Test Number: 6904767

Revision Date: 05/18/2020

TEST NAME	METHODOLOGY	LOINC CODE
Arbovirus Antibodies, IgG and IgM, Serum	Semi-Quant IFA/Semi-Quant ELISA	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL (2 mL)	Serum	Clot Activator SST	Frozen
Instructions	<p>Specimen Type: Gel-barrier tube</p> <p>Specimen Storage: FREEZE</p> <p>Specimen Collection: Serum must be separated from cells within 45 minutes of venipuncture. Send serum in a plastic transfer tube.</p> <p>Special Instructions: Testing referred to BioAgilytix Diagnostics Lab E=#-</p> <p>Specimen Stability: Ambient: 48 hours, Refrigerated : 30 days, Frozen: 30 days</p>			

GENERAL INFORMATION	
Expected TAT	6-9 Days
Performing Labcorp Test Code	820466
Notes	Labcorp Test Code: 820466
CPT Code(s)	86651x2 86652x2 86653x2 86654x2
Lab Section	Reference Lab