

Cadmium, Urine

Order Name: Cadmium Ur

Test Number: 5613555

Revision Date: 12/10/2022

TEST NAME	METHODOLOGY	LOINC CODE
Cadmium, Urine		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	8 mL (1 mL)	Urine, Random	Sterile Urine container	Room Temperature
Instructions	<p>Specimen Type: Urine Bottle</p> <p>Specimen Storage: Maintain specimen at room temperature.</p> <p>Specimen Collection: Sampling time is not critical for industrial exposure monitoring. Metals with timing "not critical" have very long elimination half-lives and accumulate in the body over years, some for a lifetime. After a couple of weeks of exposure, specimens can be collected at any time.</p> <p>Special Instructions: This test may exhibit interference when sample is collected from a person who is consuming a supplement with a high dose of biotin (also termed as vitamin B7 or B8, vitamin H, or coenzyme R). It is recommended to ask all patients who may be indicated for</p> <p>Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	2 - 4 days
Performing Labcorp Test Code	072249
Notes	Labcorp Test Code: 072249
CPT Code(s)	82300, 82570
Lab Section	Reference Lab