

C-Reactive Protein Inflammation

Order Name: **CRPINFLA**
Revision Date: 11/13/2019

| TEST NAME | | METHODOLOGY | | LOINC CODE |
|---------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|--------------------|-----------------------|
| C-Reactive Protein Inflammation | | | | |
| | | | | |
| SPECIMEN REQUIREMENTS | | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1.0 mL (0.5 mL) | Serum | Serum Gel Tube | Refrigerated |
| Alternate 1 | 1.0 mL (0.5 mL) | Plasma | Lithium Heparin | Refrigerated |
| Instructions | Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered. | | | |
| Reference Range | See Interpretive Report | | | |
| Methodology | Immunoturbidimetric | | | |

| GENERAL INFORMATION | |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Testing Schedule | Monday through Sunday |
| Expected TAT | 3 Hours |
| Stat TAT | 45 Minutes |
| Notes | Stability: Unspun - 4 Hours Room Temp - 11 Days Refrigerated - 2 Months Frozen (Serum Only) - 3 Months Do Not Freeze Plasma |
| Lab Section | NRLS-Core Chemistry |