Labcorp Oklahoma, Inc. Test Directory

Dihydrotestosterone

Order Name: DHTS

Revision Date: 10/10/2019

TEST NAME	METHODOLOGY	LOINC CODE
Dihydrotestosterone		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.0 mL (0.6 mL)	Serum	Serum Gel Tube	Refrigerated
Alternate 1	1.0 mL (0.6 mL)	Serum	Red Top	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	See interpretive report.			
Methodology	Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)			

GENERAL INFORMATION		
Testing Schedule	Monday, Wednesday, Thursday, Friday	
Expected TAT	3 Days	
Notes	Specimen Stability: Refrigerated (preferred) - 7 Days Frozen - 90 Days	
CPT Code(s)	82642 G0480 (if appropriate)	
Lab Section	NRLS-Mayo Medical Laboratories	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2024