

Lead, Capillary, with Demographics

Order Name: **PBDC**
Revision Date: 10/22/2019

TEST NAME		METHODOLOGY		LOINC CODE
Lead, Capillary, with Demographics				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.4 mL (0.1 mL)	Whole Blood	Lavender Top K2 EDTA	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			
Methodology	Inductively Coupled Plasma-Mass Spectrometry (ICP-MS)			

GENERAL INFORMATION	
Testing Schedule	Monday-Friday 5 p.m.; Saturday 2 p.m.
Expected TAT	2 Days
Notes	Specimen Stability: Refrigerated (preferred) - 28 Days Ambient - 28 Days Frozen - 28 Days
CPT Code(s)	83655
Lab Section	NRLS-Mayo Medical Laboratories