

Testosterone, Total and Bioavailable

Order Name: **TTBS**
Revision Date: 10/24/2019

TEST NAME		METHODOLOGY		LOINC CODE
Testosterone, Total and Bioavailable				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.0 mL (0.6 mL)	Serum	Red Top	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			
Methodology	TTST (Total Testosterone): Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)			
	BATS (Bioavailable Testosterone): Differential Precipitation			

GENERAL INFORMATION	
Testing Schedule	Monday-Friday
Expected TAT	2 Days
Notes	Specimen Stability: Refrigerated (preferred) - 14 Days Frozen - 60 Days
CPT Code(s)	84403 84410
Lab Section	NRLS-Mayo Medical Laboratories