

Testosterone, Total by Mass Spectrometry

Order Name: **TTST**  
Revision Date: 10/24/2019

TEST NAME		METHODOLOGY		LOINC CODE
Testosterone, Total by Mass Spectrometry				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.0 mL (0.215 mL)	Serum	Red Top	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			
Methodology	Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS) Portions of this test are covered by patent(s) held by Quest Diagnostics			

GENERAL INFORMATION	
Testing Schedule	Monday through Saturday; Continuous until 2 p.m.
Expected TAT	3 Days
Notes	Specimen Stability: Refrigerated (preferred) - 14 Days Frozen - 60 Days
CPT Code(s)	84403
Lab Section	NRLS-Mayo Medical Laboratories