Labcorp Oklahoma, Inc. **Test Directory**

Vasoactive Intestinal Polypeptide

Order Name: VIP

Revision Date: 10/24/2019

TEST NAME	METHODOLOGY	LOINC CODE
Vascactive Intestinal Polynentide		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.0 mL (0.55 mL)	Plasma	Lavender Top K2 EDTA	Frozen
Instructions	Patient must be fasting (8 hours) Spin specimen immediately and freeze. Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provi	ded.		
Methodology	Radioimmunoassay (RIA)			

GENERAL INFORMATION		
Testing Schedule	Monday, Tuesday, Wednesday; 2 p.m.	
Expected TAT	4 Days	
Notes	Specimen Stability: Frozen - 90 Days	
CPT Code(s)	84586	
Lab Section	NRLS-Mayo Medical Laboratories	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2024