Labcorp Oklahoma, Inc. Test Directory

Ehrlichia chaffeensis (HME) Antibody, IgG

Order Name: EHRC

Revision Date: 12/05/2019

TEST NAME	METHODOLOGY	LOINC CODE	
Ehrlichia chaffeensis (HME) Antibody, IgG			

SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	0.5 mL (0.15 mL)	Serum	Serum Gel Tube	Refrigerated	
Alternate 1	0.5 mL (0.15 mL)	Serum	Red Top	Refrigerated	
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.				
Reference Range	<1:64 Reference values apply to all age	es.			
Methodology	Immunofluorescence Assay (IFA)				

GENERAL INFORMATION		
Testing Schedule	Monday through Friday; 9 a.m.	
Expected TAT	1-3 Days	
Notes	Specimen Stability: Refrigerated(preferred) - 14 Days Frozen - 14 Days	
CPT Code(s)	86666	
Lab Section	NRLS-Mayo Medical Laboratories	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2024