

Parvovirus B19 Antibodies, IgG and IgM

Order Name: **PARVS**
Revision Date: 12/09/2019

TEST NAME	METHODOLOGY	LOINC CODE
Parvovirus B19 Antibodies, IgG and IgM		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.5 mL (0.5 mL)	Serum	Serum Gel Tube	Refrigerated
Alternate 1	0.5 mL (0.5 mL)	Serum	Red Top	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	IgG: Negative IgM: Negative			
Methodology	Enzyme Immunoassay (EIA)			

GENERAL INFORMATION	
Testing Schedule	Monday through Friday; 11 a.m.
Expected TAT	1-3 days
Notes	Specimen Stability: Refrigerated (preferred) - 14 Days Frozen - 14 Days
CPT Code(s)	86747 x 2
Lab Section	NRLS-Mayo Medical Laboratories