T4 (Thyroxine), Total Only

Order Name: **T4** Revision Date: **12/16/2019**

TEST NAME		MET	HODOLOGY	LOINC CODE
T4 (Thyroxine), Total Only				
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.0 mL (0.5 mL)	Serum	Serum Gel Tube	Refrigerated
Alternate 1	1.0 mL (0.5 mL)	Serum	Red Top	Refrigerated
Instructions	Patient Preparation: For 12 hours before this test do not take multivitamins or dietary supplements containing biotin (vitamin B7), which is commonly found in hair, skin, and nail supplements and multivitamins.			
	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			
Methodology	Electrochemiluminescence Immunoassay			
GENERAL INFORMATION				
Testing Schedule	Monday through Friday 5 a.m 12 a.m.; Saturday 6 a.m 6 p.m.			
Expected TAT	1-3 Days			
Notes	Specimen Stability: Refrigerated (preferred)- 7 Days Frozen - 30 Days			
CPT Code(s)	84436			
Lab Section	NRLS-Mayo Medical Laboratories			