

Trichomonas vaginalis, Miscellaneous Sites, by Nucleic Acid Amplification

Order Name: **MTRNA**
Revision Date: 05/21/2020

TEST NAME	METHODOLOGY	LOINC CODE
Trichomonas vaginalis, Miscellaneous Sites, by Nucleic Acid Amplification		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	15-20 mL (15-20 mL)	Urine, Random	APTIMA Urine Transport Tube	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	Negative			
	Interpretation:A positive result is considered indicative of current or recent Trichomonas vaginalis infection (trichomoniasis).			
Methodology	Transcription Mediated Amplification			

GENERAL INFORMATION	
Testing Schedule	Monday through Saturday; Varies
Expected TAT	1-3 Days
Notes	Specimen Stability: Refrigerated: 30 Days Frozen- 30 Days
CPT Code(s)	87661
Lab Section	NRLS-Mayo Medical Laboratories