

## Parainfluenza 1,2,3 by Complement Fixation

Order Name: **PARAINFLU**  
Test Number: 5504941  
Revision Date: 12/15/2022

| TEST NAME                                  | METHODOLOGY              | LOINC CODE |
|--|--------------------------|------------|
| Parainfluenza 1,2,3 by Complement Fixation | Complement fixation (CF) | 1668-3     |

### SPECIMEN REQUIREMENTS

|              |   |               |                           |                       |
|--------------|---|---------------|---------------------------|-----------------------|
| Specimen     | Specimen Volume (min)   | Specimen Type | Specimen Container        | Transport Environment |
| Preferred    | <b>2 mL</b>   | <b>Serum</b>  | <b>Clot Activator SST</b> | <b>Refrigerated</b>   |
| Instructions | <b>Specimen Type:</b> Red-top tube or gel-barrier tube<br><b>Specimen Storage:</b> Refrigerate<br><b>Specimen Collection:</b> Not Available<br><b>Specimen Stability:</b> Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available |               |                           |                       |

### GENERAL INFORMATION

|                              |                    |
|------------------------------|--------------------|
| Expected TAT                 | 1 - 5 days 2 days. |
| Performing Labcorp Test Code | 096214             |
| CPT Code(s)                  | 86710x3            |
| Lab Section                  | Reference Lab      |