

## HemePath Consult (w/CBC)

Order Name: **HemePath Consult**  
Test Number: 2904600  
Revision Date: 11/16/2023

TEST NAME	METHODOLOGY	LOINC CODE
<b>Complete Blood Count (CBC) with Automated Differential</b>	Flow cytometry	See Individual Assays
<b>Reticulocyte (Retic) Count</b>	Flow cytometry	
<b>Immature Platelet Fraction</b>	Flow cytometry	71693-6
<b>HemePath Consult order</b>		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
<b>Preferred</b>	<b>5 mL (1 mL)</b>	<b>Whole Blood</b>	<b>EDTA (Lavender Top)</b>	<b>Refrigerated</b>
<b>Alternate 1</b>	<b>1 mL (0.5 mL)</b>	<b>Whole Blood</b>	<b>EDTA (Lavender) Microtainer/Bullet</b>	<b>Refrigerated</b>
<b>Instructions</b>	<p><b>For best results: Room temperature specimens should be tested within 24hrs, otherwise send Refrigerated.</b> Refrigerated specimens can be tested up to 48hrs. Specimens received after 24hrs will not receive a manual differential. Slides must be made within 24hr of collection. <b>Specimens received greater than 48hrs old will be canceled.</b></p> <p><i>Note: This test cannot be added on to an A1c that has already been performed.</i></p>			

GENERAL INFORMATION	
<b>Testing Schedule</b>	Daily
<b>Expected TAT</b>	1-2 Days
<b>Clinical Use</b>	Pathology Report will be issued on this on Order/Result Code: <b>[8077777] Hematology Other Report</b> , at an additional charge to the patient.
<b>Notes</b>	Provide patient history as available. Testing includes CBC with automated differential, IPF, Retic Count and pathologist interpretation. If this testing is performed at your laboratory please send these results with the smears and the lavender tube. If the question is anemia, consider ordering an Anemia Analyzer with the algorithmic reflex ordering of the appropriate chemistry tests.
<b>CPT Code(s)</b>	85025, 85045, 85055, Additional professional code determined with clinical findings
<b>Lab Section</b>	Hematology