

## Haptoglobin

Order Name: **HAPTOGLOBN**

Test Number: 5000700

Revision Date: 05/20/2021

TEST NAME	METHODOLOGY	LOINC CODE
Haptoglobin	Immunoturbidimetry	4542-7

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Alternate 1	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Instructions	Stability: Ambient 4 hours, Refrigerated 7 days, Frozen 2 weeks			

### GENERAL INFORMATION

Testing Schedule	Mon - Fri
Expected TAT	1-3 days
Clinical Use	A low serum haptoglobin concentration is not specific for hemolysis. Evaluation and diagnosis of hemolytic anemia.
CPT Code(s)	83010
Lab Section	Chemistry