

Immunoglobulin IgG4 Quantitative

Order Name: **IGG4**
Test Number: 5580200
Revision Date: 10/23/2017

| TEST NAME | METHODOLOGY | LOINC CODE |
|----------------------------------|---------------|------------|
| Immunoglobulin IgG4 Quantitative | Turbidometric | 2469-5 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|--|---------------|--------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1mL | Serum | Clot Activator SST | Refrigerated |
| Instructions | Allow specimen to clot completely at room temperature. Separate serum or plasma from cells ASAP or within 2 hours of collection. Stability After separation from cells: Ambient 4hours, Refrigerated 7 days, Frozen 1 months (avoid repeated freeze/thaw cycles). | | | |

| GENERAL INFORMATION | |
|---------------------|-----------------------|
| Testing Schedule | Mon-Sat |
| Expected TAT | 2-8 Days |
| CPT Code(s) | 82787 |
| Lab Section | Immunology - Serology |