

## CA 125, Serous Fluid

Order Name: **SRS CA125**

Test Number: 2015430

Revision Date: 06/14/2018

TEST NAME	METHODOLOGY	LOINC CODE
CA 125, Serous Fluid	Chemiluminescence Assay	11210-2

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serous Fluid	Sterile Screwtop Container	Refrigerated
Instructions	Venous blood is often drawn simultaneously. Note fluid type on requisition and container. Stability: Ambient 1 day. Refrigerated 7 days. Frozen > 7 days.			

### GENERAL INFORMATION

Testing Schedule	Mon - Fri
Expected TAT	1-3 days
Clinical Use	Useful in the follow-up management of patients undergoing cancer therapy, especially for ovarian carcinoma.
CPT Code(s)	86304
Lab Section	Chemistry