

## Testosterone Total, Adult Male

Order Name: **TESTOS TO**

Test Number: 3602650

Revision Date: 06/15/2018

TEST NAME	METHODOLOGY	LOINC CODE
Testosterone Total, Adult Male	Chemiluminescence Assay	2986-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Alternate 1	1 mL (0.5)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	<p><b>Not Recommended for Females and Children 13 years and under</b></p> <p>Preferred to separate serum from cells ASAP or within 2 hours of collection. Stability: Ambient 8 hrs. Refrigerated 7 days. Frozen &gt; 7 days</p>			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	1-2 days
Clinical Use	Useful for evaluation of males with erectile dysfunction, gynecomastia, osteoporosis, infertility, delayed or precocious puberty and monitoring replacement therapy. Evaluation of women with hirsutism, virilization and oligomenorrhea and infants with ambiguous genitalia and/or virilizing syndromes.
Notes	<p><b>For Females and Children 13 years and under use the following test code:</b> Testosterone Total by LCMS [3602651]</p>
CPT Code(s)	84403
Lab Section	Chemistry