

## Gastrin

Order Name: **GASTRIN**  
Test Number: 3601300  
Revision Date: 12/10/2022

TEST NAME	METHODOLOGY	LOINC CODE
Gastrin	Immunochemiluminometric (ICMA)	2333-3

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.5 mL (0.3 mL)	Serum	Clot Activator SST	Frozen

**Instructions**

**Notes:** 0.3 mL (Note: This volume Does NOT allow for repeat testing.)  
**Specimen Type:** Red-top tube or gel-barrier tube  
**Container Detail:** LC PP  
**Specimen Storage:** FREEZE immediately.  
**Specimen Collection:** Separate serum from cells. Transfer the serum into a LabCorp PP transpak frozen purple tube with screw cap (LabCorp ID 49482). Freeze immediately and maintain frozen until tested. To avoid delays in turnaround time when requesting multiple test on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.  
**Special Instructions:** This test may exhibit interference when sample is collected from a person who is consuming a supplement with a high dose of biotin (also termed as vitamin B7 or B8, vitamin H, or coenzyme R). It is recommended to ask all patients who may be indicated for  
**Specimen Stability:** Ambient: Not Available, Refrigerated : Not Available, Frozen: 14 days

### GENERAL INFORMATION

Expected TAT	2-4 Days
Clinical Use	Aids in diagnosis of carcinoid and gastrinoma tumors.
Performing Labcorp Test Code	004390
Notes	Labcorp Test Code: 004390
CPT Code(s)	82941
Lab Section	Reference Lab