

Electrolytes, Feces

Order Name: **ELECT FEC**

Test Number: 3605725

Revision Date: 12/12/2022

| TEST NAME | METHODOLOGY | LOINC CODE |
|-----------------|-------------------------|------------|
| Potassium Stool | Ion-Selective Electrode | 15202-5 |
| Sodium Stool | Ion-Selective Electrode | 15207-4 |
| Chloride Stool | Ion-Selective Electrode | 15158-9 |

SPECIMEN REQUIREMENTS

| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
|--------------|---|------------------|---|-----------------------|
| Preferred | 10gm (1gm) | Timed Stool 24hr | Stool Specimen Container, Large for 24 or 72 hour Collections (White) | Frozen |
| Alternate 1 | 10gm (1gm) | Stool, Random | Stool specimen container | Frozen |
| Instructions | <p>Specimen Type: Sterile container Specimen Storage: Refrigerated Specimen Collection: Stool must be liquid. Do not use saline or water to liquify sample. Special Instructions: Testing referred to ARUP TC 60185 Specimen Stability: Ambient: Unacceptable, Refrigerated : 1 week, Frozen: 1 month</p> | | | |

GENERAL INFORMATION

| | |
|------------------------------|---------------------------|
| Expected TAT | 4 - 8 days |
| Performing Labcorp Test Code | 827678 |
| Notes | Labcorp Test Code: 827678 |
| CPT Code(s) | 82438; 84999; 84302 |
| Lab Section | Reference Lab |