

Lorazepam (Ativan)

Order Name: **LORAZEPAM**

Test Number: 3611975

Revision Date: 12/20/2022

TEST NAME	METHODOLOGY	LOINC CODE
Lorazepam (Ativan)	Quantitative Liquid Chromatography/Tandem Mass Spectrometry	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3 mL (0.6 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen Collection: Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube. Special Instructions: Testing referred to MedTox Laboratories Inc MX#64			

GENERAL INFORMATION

Expected TAT	7-10 Days
Performing Labcorp Test Code	808465
Notes	Labcorp Test Code: 808465
CPT Code(s)	80299
Lab Section	Reference Lab