## Labcorp Oklahoma, Inc. Test Directory

## Actin (F-Actin) Smooth Muscle Antibody IgG

Order Name: ACTIN AB Test Number: 5700200 Revision Date: 12/12/2022

TEST NAME			METHODOLOGY		
Actin (F-Actin) Smooth Muscle Antibody IgG		Enzyme-Linked Immunosorbent Assay (ELISA	44706-0		
SPECIMEN REQUIREM	ENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	1 mL (0.3)	Serum	Clot Activator SST	Room Temperature	
Instructions	Specimen Type: Red-top tube Specimen Storage: Room tem Specimen Collection: Not Ava Special Instructions: This ass	becimen Type: Red-top tube or gel-barrier tube   becimen Storage: Room temperature   becimen Collection: Not Available   becimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days			
GENERAL INFORMATIO	ON				
Expected TAT	2-3 Days				
Clinical Use	In viral hepatitis the tit to 1:40.	In viral hepatitis the titers are generally less than 1:80 and are transient. The titers in primary biliary cirrhosis are also low, ranging from 1:20 to 1:40.			
Performing Labcorp Tes	t 006643				

Code			
Notes	Labcorp Test Code: 006643		
CPT Code(s)	86015		
Lab Section	Reference Lab		