

Helicobacter pylori Antigen Detection Stool

Order Name: **HPYLORI AG**

Test Number: 3502325

Revision Date: 09/21/2022

TEST NAME	METHODOLOGY	LOINC CODE
Helicobacter pylori Antigen Detection Stool	Immunoassay (IA)	80373-4

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5g (1)	Stool, Random	Stool specimen container	Refrigerated

Instructions

Patient Preparation: To avoid false negative results, no antibiotics, proton pump inhibitors, or a bismuth preparation should be administered to patient for 14 days prior to specimen collection.

Transfer 5g(1g) stool to an unpreserved stool transport vial and Freeze.

Unacceptable Conditions: Gastric specimens, swabs or tissue. Specimens in media or preservatives.

Stability: Ambient: 2 hours; Refrigerated: 72 hours; Frozen: 1 month.

GENERAL INFORMATION

Testing Schedule	Sun-Sat
Expected TAT	2 Days
CPT Code(s)	87338
Lab Section	Microbiology