

Eye Culture and Stain

Order Name: **C EYE RTS**

Test Number: 6002006

Revision Date: 09/21/2017

| TEST NAME | METHODOLOGY | LOINC CODE |
|-----------------------|-------------|------------|
| Eye Culture and Stain | Culture | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|------------------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | See Instructions | Swab | Copan eSwab - White (Regular size) | Room Temperature |
| Alternate 1 | See Instructions | Swab | Copan eSwab - Green (Mini-tip) | Room Temperature |
| Alternate 2 | See Instructions | Swab | PNP Swab (Green Cap) | Room Temperature |
| Instructions | Swab specimen should remain at room temperature and be setup within 2 days of collection. Dry swabs, swabs in non-sterile transport devices, viral media, or refrigerated and/or frozen swabs are unacceptable. | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Testing Schedule | Daily |
| Expected TAT | 3 Days |
| Clinical Use | Eye cultures help determine bacterial and fungal infections in and around the eye. Swab specimens should be carefully collected to maintain sterility. |
| CPT Code(s) | 87070 |
| Lab Section | Microbiology |